

Do hypoxia/normoxia culturing conditions change the neuroregulatory profile of Wharton Jelly mesenchymal stem cell secretome?

Fábio G. Teixeira, Krishna M. Panchalingam, Sandra Isabel Anjo, Bruno Manadas, Ricardo Pereira, Nuno Sousa, António J. Salgado, and Leo A. Behie

Presented by:


Suzana Tufegdzcic

Stem Cells


- undifferentiated cells that are able to differentiate into specialized cells
- Classical abilities: self-renewel and potency
- Two types of stem cells:
 - **embryonic stem cells** → isolated from blastocysts (4-5 days post fertilization)
 - **adult stem cells** → in various tissues → "repair system of the body"
 - Sources of autologous adult stem cells: bone marrow, adipose tissue, blood, umbilical cord
- Treatment:
 - Bone marrow transplant → for patients with cancers of blood or bone marrow (multiple myeloma or leukemia)
 - Problems: immunosuppression, no specific cell type → pluripotency, some stem cells form tumors

Stem cells

Research



Wharton Jelly Mesenchymal Stem Cells



Cross-section of an umbilical cord


Wharton Jelly → gelatinous substance within the umbilical cord

Mesenchymal Stem Cells (MSC) → able to differentiate into bone, nerve, cornea, heart, fat and cartilage cells

→ sources of MSC: umbilical cord,


bone marrow
adipose
tissue
tooth pulp

Reaction of cells in hypoxic conditions:



- most important parameter in healthy cells are a high amount of ATP
- Hypoxia → Mitochondria produces less and less ATP → failure of Na⁺/K⁺ ATPase → membrane depolarization → Ca²⁺ influx → activation of calcium dependent-phospholipases and proteases → cell swelling, hydrolysis of cellular components

Effect of hypoxia on gene expression:



- Hypoxia inhibits the oxygen sensor PHD and FIH → HIF- α is stable
→ translocates to the nucleus → dimerizes with HIF- β → binds to DNA at a HRE (hormone response element) → activates the transcription of genes

PHD = Prolyl-4-hydroxylase

FIH = Factor inhibiting HIF

HIF- α = Hypoxia-inducible factor 1- α

VEGF = vascular endothelial growth factor

Inos2 = inducible nitric oxide synthase

ho1 = haemoxgenase-1

epo = erythropoetin

Hypoxic Conditions



Introduction

- human Mesenchymal Stem Cells (hMSC) → important role in regenerative medicine
- The problem → low amount of cells after isolation → difficult to apply in clinical conditions
- hMSC → normal in static culture and in medium with 21% oxygen tension
- BUT hMSC are in the body at much lower oxygen tension (bone marrow 1-7%, adipose tissue 10-15%, umbilical cord 1,5-8%)
- The study from Rhijn et al. (2013) showed that hypoxic conditioning of hMSC improve the regenerative potential and so preserve their immunosuppressive ability.
- The study from Tsai et al. (2011) used 1 % oxygen that showed the reduction of their senescence and augmentation of their proliferation

Introduction

- angiogenic potential of MSCs is affected by modulating the oxygen concentration → a higher level of VEGF, bFGF and HGF

(**VEGF**=vascular endothelial growth factor, **bFGF** = beta-fibroblast growth factor, **HGF** =hepatocyte growth factor)

- cell death at long exposure to hypoxic conditions (Volkmer et al., 2008)
- Higher oxygen tension might be toxic → oxidative stress caused by production of ROS (=reactive oxygen species)

Methods

- DASGIP Parallel Bioreactor System:
 - growth of human Wharton Jelly Mesenchymal Stem Cells
 - contains Teflon 4-paddle impellers
 - siliconized with Sigmacote and autoclaved of vessels
 - Hold at 37 °C with a heating jacket
 - 100% dissolved oxygen for normoxic conditions
 - 21% dissolved oxygen for the hypoxic conditions
 - pH =7,4
 - Agitation of 52 rpm using a magnetic stir plate under the bioreactor



Methods

microcarriers

- utilizing Cytodex 3 microcarriers
- 1,0 g (per condition) of microcarriers hydrated in 50 ml 1x phosphate-buffered saline (PBS) in 2 x 125ml Erlenmeyer flasks at room temperature for one night
- 2-3 drops of Tween 80 → for breaking surface tension and guarantee proper wetting and sedimentation
- 3 x times washing of microcarriers with 1xPBS and autoclaved
- incubated with fetal bovine serum (FBS) for 6 h at 37 °C for coating the microcarriers with serum-attachment factors

Methods

Preparation of microcarriers:

- after 6h, microcarriers washed with serum-free medium (PPRF-msc6)
- cocultured in 275 ml of serum-free medium for 4 h at controlled culture conditions

Methods

- hWJ-MSCs:
 - expanded in Nunc T-flask with the gelatine solution (including „ Typ B bovine gelatine“)
 - In static culture for two passages
 - harvested with trypsin-EDTA and inoculate in bioreactors
 - Volume of the bioreactors →325ml for the first 24 h and then to 500ml increased
 - Bioreactors removed from DASGIP system and placed for 10 min a biosafety cabinet
 - Cells were cultured on the microcarriers for 72 h!

Methods

- the supernatant was removed from the bioreactors
- microcarriers washed with Neurobasal[®]-A medium
- 500ml of Neurobasal[®]-A medium with 1 % kanamycine were added in the bioreactors
- again for 24 h in the DASGIP control system
- then removed from the system and for 10 min in the biosafety cabinet
- supernatant was harvested →centrifuged at 300g for 10 min to take off the cell debris

Methods

Human telencephalon neural precursor cell (hNPCs)

- hNPCs
 - isolated from telecephalon area of a 10-week fetus and cryopreserved
 - melted at 37°C and placed into a Nunc T-25 flask including 5 ml of serum-free medium (Δ PPRF-h2)
 - after 2 day \rightarrow cell collected and mechanically dissociated into single cell suspension
 - every 4 days 40 % of the medium was replaced by fresh growth medium
 - after 14-20 days hNPCs were passaged and plated into precoated laminin 24-well plates including the hWJ-MSc condition medium from normoxic and hypoxic condition
 - Neurobasal[®]-A was the control group
 - the plates \rightarrow placed in a incubator with 37°C, 5% CO₂, 95% air and 90% relative humidity

Methods

Immunostaining

- hNPCs:
 - fixed in 4% paraformaldehyde for 15 min
 - permeabilized by incubation with 0,1% Triton X-100 in 1xPhosphate buffered saline (PBS)
 - washed 3 times in 1xPBS
 - blocked with 10 % fetal calf serum in 1xPBS
 - 1 h incubation with primary antibodies
 - rabbit anti-doublecortin (DCX) → detect immature neurons
 - mouse antiratmicrotubule associated protein-2 (MAP-2) → detect mature neurons

Methods

Immunostaining

- hNPCs:
 - 3 times washed in 1xPBS
 - incubated with the secondary antibodies: Alexa Fluor 488 goat anti-rabbit and Alexa Fluor 594 goat anti-mouse immunoglobulin G for 1 h at 37°C
 - 10 min with 4-6-diamidino-2-phenylindole-dihydrochloride
 - stained and then observed under an Olympus BX-61 Fluorescence Microscope
 - normalizing the data with presenting the results in percentage of cells → calculated by counting the cells positive for MAP-2/DCX markers and dividing it by total number of cells/field

Methods

Sample preparation and quantification

- Vivaspin 20 sample concentrator → ultracentrifugation at 3000g for 45 min
- Condition medium → precipitated with Trichochloroacetic – acetone
- Samples were evaluated with 2D-Quant Kit
- Samples were digested with trypsin at 37°C
- peptides were dried by rotary evaporation
- Triple TOF 5600System → analyzed the samples (mass spectrometry)

Results / Discussion




Fig.1 (A,B) → hWJ-MSCs could attach to and grow on Cytodex 3 microcarriers

Results / Discussion







Fig.2 → expansion of immature and mature neurons in hypoxic and normoxic conditions

Results / Discussion



Venn diagram:
 → more proteins in in hWJ-MSC
 CM from hypoxic conditions



Red → low expression
Green → high expression

Fig.3 → WJ-MSC CM proteomic analysis by mass spectrometry

Results / Discussion




Fig. 4 → both hWJ-MSC secretomes were analyzed for specific proteins with neuroregulatory potential

Results / Discussion

Normoxic conditions

- Vironectin
- Cadherin-2
- MRP-1

Hypoxic conditions

- PEDF
- IGF-2
- MIF
- Hsp 70
- semaphorin-7A
- moesin

Results/ Discussion

- hypothesized normoxic conditions could cause a reduction in therapeutic potential of human Mesenchymal Stem Cells
→in vivo physiologic oxygen concentration is lower than 21 % oxygen tension
- this study showed in both conditions that hWJ-MSCs can be cultivated and induce differentiation of hNPCs


Discussion

- Pro:
 - reproducibility → very detailed methods
 - illustrations understandable, except of Fig. 1 → C, D
- Cons:
 - lack of definition of mature and immature neurons
 - Glia cells ?

Discussion:

- Other studies:
 - from Rhijn et al. (2013) → “Effects of hypoxia on immunomodulatory properties of adipose tissue-derived mesenchymal stem cells”:
 - Effect of adipose tissue-derived mesenchymal stem cells (ASC) in hypoxic (1% oxygen) culture conditions → no immunophenotypical changes → no effect of proliferation of ASC in 1% and 20% oxygen!
 - from Cicione et al. (2013) → “Effects of Severe Hypoxia on Bone Marrow Mesenchymal Stem Cells Differentiation Potential”:
 - absence of differentiation in severe hypoxic conditions compared to normoxic condition → severe culture condition

Any questions?



Thank you for your attention!

